

## Anti-TET2 Antibody

Catalog Number: A00191-1

### About TET2

TET2, a member of the ten-eleven-translocation (TET) family of genes, is a methylcytosine dioxygenase that catalyzes the conversion of methylcytosine to 5-hydroxymethylcytosine. It is a candidate tumor suppressor gene reported to be mutated in approximately 14% of patients with JAK2V617F-positive myeloproliferative neoplasms (1), and can be mutated in other hematopoietic disorders such as myelodysplastic syndromes, acute myeloid leukemia, and chronic myelomonocytic leukemia (2). Analysis of the TET2 and JAK2 mutations in these neoplasms suggests that mutations in TET2 do not represent a predisposition for acquiring mutations in JAK2.

### Overview

Product Name	Anti-TET2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-TET2 Antibody (Catalog # A00191-1). Tested in ELISA, WB applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	TET2 Antibody is supplied in PBS containing 0.02% sodium azide.
Storage Instructions	TET2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Host	Rabbit
Uniprot ID	Q6N021

### Technical Details


Immunogen	TET2 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human TET2. The immunogen is located within amino acids 1550 - 1600 of TET2.
Predicted Reactive Species	Mouse
Cross Reactivity	TET2 antibody is human and mouse reactive. At least two isoforms of TET2 are known to exist; this antibody will detect the larger isoforms. This antibody is predicted to not cross-react with TET1 and TET3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/mL
Purification	TET2 antibody is affinity chromatography purified via peptide column.

Suggested Dilutions

TET2 antibody can be used for detection of TET2 by Western blot at 1 - 2 ug/ml.  
Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested. Optimal dilutions for each application should be determined by the researcher.

## Anti-TET2 Antibody (A00191-1) Images

---



Western blot analysis of TET2 in 3T3 cell lysate with TET2 antibody at 1 ug/ml.

---

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TET2 Antibody

For Research Use Only. Not for use in diagnostic procedures.