

Anti-FGFR3 Antibody

Catalog Number: A00200-1


Overview

Product Name	Anti-FGFR3 Antibody
Reactive Species	Human
Description	Boster Bio Anti-FGFR3 Antibody catalog # A00200-1. Tested in WB,ICC/IF,IHC,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P22607


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-FGFR3 Antibody (A00200-1) Images



Western blot analysis of FGFR3 on HeLa cells lysates using anti-FGFR3 antibody at 1/500 dilution.



ICC staining FGFR3 in HeLa cells (red). The nuclear counterstain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilized with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 33626078, Nannapaneni S,Griffith CC,Magliocca KR,Chen W,Lyu X,Chen Z,Wang D,Wang X,Shin DM,Chen ZG,Saba NF.Co-expression of fibroblast growth factor receptor 3 with mutant p53, and its association with worse outcome in oropharyngeal squamous cell carcinoma.PLoS One

Visit bosterbio.com/anti-fgfr3-antibody-a00200-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGFR3 Antibody

For Research Use Only. Not for use in diagnostic procedures.