

## Anti-Rho A (R182) Antibody

Catalog Number: A00207

### Overview


Product Name	Anti-Rho A (R182) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Rho A (R182) Antibody catalog # A00207. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P61586

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 150-200 of Human Rho A.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

## Anti-Rho A (R182) Antibody (A00207) Images

---



Immunohistochemistry (IHC) analyzes of Rho A (R182) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

## 1 Publications Citing This Product

---

1. PubMed ID: 10.12659/MSM.907277, Rho Kinase Type 1 (ROCK1) Promotes Lipopolysaccharide-induced Inflammation in Corneal Epithelial Cells by Activating Toll-Like Receptor 4 (TLR4)-Mediated Signaling

Visit [bosterbio.com/anti-rho-a-r182-antibody-a00207-boster.html](https://bosterbio.com/anti-rho-a-r182-antibody-a00207-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Rho A (R182) Antibody

For Research Use Only. Not for use in diagnostic procedures.