

Anti-IL-2R Alpha (R264) IL2RA Antibody

Catalog Number: A00214


Overview

Product Name	Anti-IL-2R Alpha (R264) IL2RA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IL-2R Alpha (R264) IL2RA Antibody catalog # A00214. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01589

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 231-280 of Human IL-2Ralpha.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IF: 1:50-1:200

Anti-IL-2R Alpha (R264) IL2RA Antibody (A00214) Images



Western blot (WB) analysis of IL-2Ralpha (R264) polyclonal antibody at 1:500 dilution Lane1:PC3 whole cell lysate(40ug) Lane2:A549 whole cell lysate(40ug) Lane3:PMVEC whole cell lysate(40ug) Lane4:BV2 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-2R Alpha (R264) IL2RA Antibody

For Research Use Only. Not for use in diagnostic procedures.