

## Anti-Axl Antibody

Catalog Number: A00226-3


### Overview

Product Name	Anti-Axl Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Axl Antibody catalog # A00226-3. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P30530

### Technical Details

Immunogen	Synthesized peptide derived from Axl . at AA range: 630-710
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-Axl Antibody (A00226-3) Images



Western blotting validation for Anti-Axl Antibody A00226-3  
Western Blot (WB) analysis of 293T using Axl antibody.  
Electrophoresis was performed on a SDS-PAGE gel. To  
determine SDS-PAGE gel concentration

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-Axl Antibody

For Research Use Only. Not for use in diagnostic procedures.