

Anti-IL-4 Antibody

Catalog Number: A00230-2


Overview

Product Name	Anti-IL-4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IL-4 Antibody catalog # A00230-2. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P05112

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 108-153 of Human IL-4.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-IL-4 Antibody (A00230-2) Images



Western blot (WB) analysis of IL-4 polyclonal antibody at 1:500 dilution Lane1:MCF-7 cell lysate Lane2:sp2/0 cell lysate Lane3:PC12 cell lysate

1 Publications Citing This Product

1. PubMed ID: 33683631, Zhao N, Wang T, Peng L, Li Y, Zhao Y, Yu S. Attenuation of Inflammation by DJ-1 May Be a Drug Target for Cerebral Ischemia-Reperfusion Injury. *Neurochem Res.* 2021 Mar 8. doi: 10.1007/s11064-021-03288-z. Epub ahead of print. PMID: 33683631.

Visit bosterbio.com/anti-il-4-antibody-a00230-2-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL-4 Antibody

For Research Use Only. Not for use in diagnostic procedures.