

Anti-HLA-DQA1 Antibody

Catalog Number: A00232-1

Overview

Product Name	Anti-HLA-DQA1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HLA-DQA1 Antibody catalog # A00232-1. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01909

Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 IHC: 1:50-1:200 IP: 1:10-1:50

Anti-HLA-DQA1 Antibody (A00232-1) Images




Figure 2. Immunohistochemistry validation of HLA-DQA1 using Anti-HLA-DQA1 Antibody (A00232-1).

Immunohistochemical analysis of paraffin-embedded rat lung tissue using anti-HLA-DQA1 antibody. Counter stained with hematoxylin.

For more protocol information of IHC




Figure 1. Western blotting validation for Anti-HLA-DQA1 Antibody A00232-1

Western blot analysis of HLA-DQA1 on different lysates using anti-HLA-DQA1 antibody at 1/1

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HLA-DQA1 Antibody