

Anti-FADD (N188) Antibody

Catalog Number: A00237-2


Overview

| | |
|----------------------|--|
| Product Name | Anti-FADD (N188) Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-FADD (N188) Antibody catalog # A00237-2. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat. |
| Application | IHC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q13158 |


Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic peptide, corresponding to amino acids 150-200 of Human FADD. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 IHC: 1:50-1:200 |

Anti-FADD (N188) Antibody (A00237-2) Images



Western blot (WB) analysis of FADD (N188) polyclonal antibody at 1:500 dilution Lane1:A549 whole cell lysate(40ug) Lane2:HCT116 whole cell lysate(40ug) Lane3:Panc1 whole cell lysate(40ug) Lane4:C6 whole cell lysate(40ug) Lane5:MEF whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of FADD (N188) pAb in paraffin-embedded human lung carcinoma tissue at 1:100.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FADD (N188) Antibody

For Research Use Only. Not for use in diagnostic procedures.