

Anti-Autophagy protein 5 ATG5 Antibody

Catalog Number: A00240-1


Overview

Product Name	Anti-Autophagy protein 5 ATG5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Autophagy protein 5 ATG5 Antibody catalog # A00240-1. Tested in WB,ICC/IF,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H1Y0


Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:100-1:500 IHC: 1:50-1:200

Anti-Autophagy protein 5 ATG5 Antibody (A00240-1) Images



Western blot analysis of ATG5 on different lysates using anti-ATG5 antibody at 1/1,000 dilution. Positive control: Lane 1: Raji Lane 2: HeLa Lane 3: PC12



ICC staining ATG5 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Autophagy protein 5 ATG5 Antibody

For Research Use Only. Not for use in diagnostic procedures.