

## Anti-ARK-1 AURKA Antibody

Catalog Number: A00246-2

### Overview


Product Name	Anti-ARK-1 AURKA Antibody
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-ARK-1 AURKA Antibody catalog # A00246-2. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat, Monkey.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O14965

### Technical Details


Immunogen	The antiserum was produced against synthesized peptide derived from human AurA. AA range:311-360
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:40000



## Anti-ARK-1 AURKA Antibody (A00246-2) Images



Western blotting validation for Anti-ARK-1 AURKA Antibody A00246-2 Western Blot (WB) analysis of specific cells using ARK-1 polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of AURKA using Anti-ARK-1 AURKA Antibody (A00246-2). Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4°C)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ARK-1 AURKA Antibody

For Research Use Only. Not for use in diagnostic procedures.