

Anti-HSF1 (phospho-S303) Antibody

Catalog Number: A00250S303


Overview

Product Name	Anti-HSF1 (phospho-S303) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HSF1 (phospho-S303) Antibody catalog # A00250S303. Tested in IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00613

Technical Details

Immunogen	Synthetic phosphopeptide derived from human HSF1 around the phosphorylation site of Serine 303.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 IF: 1:50-1:200

Anti-HSF1 (phospho-S303) Antibody (A00250S303) Images



Immunohistochemistry (IHC) analyzes of p-HSF1 (S303) pAb in paraffin-embedded human breast carcinoma tissue at 1:50, showing cytoplasmic and nucleus staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSF1 (phospho-S303) Antibody

For Research Use Only. Not for use in diagnostic procedures.