

Anti-HSF1 (S303) Antibody

Catalog Number: A00250S303-1


Overview

Product Name	Anti-HSF1 (S303) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-HSF1 (S303) Antibody catalog # A00250S303-1. Tested in WB applications. This antibody reacts with Human,Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q00613

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 270-315 of Human HSF1.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-HSF1 (S303) Antibody (A00250S303-1) Images



Western blot (WB) analysis of HSF1 polyclonal antibody at 1:1000 dilution Lane1:HELA whole cell lysate(41ug)
Lane2:DLD whole cell lysate(37ug) Lane3:NIH3T3 whole cell lysate(48ug) Lane4:H9C2 whole cell lysate(37ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSF1 (S303) Antibody

For Research Use Only. Not for use in diagnostic procedures.