

## Anti-Histone deacetylase 1 HDAC1 Antibody

Catalog Number: A00256-3


### Overview

Product Name	Anti-Histone deacetylase 1 HDAC1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone deacetylase 1 HDAC1 Antibody catalog # A00256-3. Tested in WB,ICC/IF,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13547


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1.33 mg/mL by UV absorbance at 280 nm
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

## Anti-Histone deacetylase 1 HDAC1 Antibody (A00256-3) Images



Western blot analysis of HDAC1 on different lysates using anti-HDAC1 antibody at 1/1,000 dilution. Positive control:  
Lane 1: Hela Lane 2: Jurkat Lane 3: K562



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using anti-HDAC1 antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone deacetylase 1 HDAC1 Antibody

For Research Use Only. Not for use in diagnostic procedures.