

Anti-NFkB p65 (Phospho-S529) RELA Antibody

Catalog Number: A00284S529-1


Overview

Product Name	Anti-NFkB p65 (Phospho-S529) RELA Antibody
Reactive Species	Human
Description	Boster Bio Anti-NFkB p65 (Phospho-S529) RELA Antibody catalog # A00284S529-1. Tested in WB,IP,Flow Cytometry applications. This antibody reacts with Human.
Application	Flow Cytometry, IP, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q04206

Technical Details

Immunogen	Synthetic peptide conjugated to Keyhole Limpet Haemocyanin (KLH)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 FC: 1:50-1:100

Anti-NFkB p65 (Phospho-S529) RELA Antibody (A00284S529-1) Images



Flow cytometric analysis of Daudi cells with Phospho-NF-kB p65(S529) antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NFkB p65 (Phospho-S529) RELA Antibody

For Research Use Only. Not for use in diagnostic procedures.