

## Anti-NFkB-p65 (T429) RELA Antibody

Catalog Number: A00284T429


### Overview

Product Name	Anti-NFkB-p65 (T429) RELA Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-NFkB-p65 (T429) RELA Antibody catalog # A00284T429. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q04206


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 400-450 of Human NFkB-p65.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-NFkB-p65 (T429) RELA Antibody (A00284T429) Images



Western blot (WB) analysis of NFkB-p65 (T429) polyclonal antibody at 1:500 dilution Lane1:CT26 whole cell lysate(40ug) Lane2:BV2 whole cell lysate(40ug) Lane3:The Kidney tissue lysate of Mouse(40ug) Lane4:NIH-3T3 whole cell lysate(40ug) Lane5:H9C2 whole cell lysate(40ug) Lane6:L02 whole cell lysate(40ug) Lane7:HepG2 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of NFkappaB-p65 (T429) pAb in paraffin-embedded human breast carcinoma tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NFkB-p65 (T429) RELA Antibody

For Research Use Only. Not for use in diagnostic procedures.