

## Anti-ErbB-4 / HER4 (phospho-Y1284) Antibody

Catalog Number: A00296Y1284


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-ErbB-4 / HER4 (phospho-Y1284) Antibody   |
| Reactive Species     | Human, Mouse, Rat   |
| Description          | Boster Bio Anti-ErbB-4 / HER4 (phospho-Y1284) Antibody catalog # A00296Y1284. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application          | WB  |
| Clonality            | Polyclonal  |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.            |
| Host                 | Rabbit  |
| Uniprot ID           | Q15303  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | Synthetic phosphopeptide derived from human ErbB-4 / HER4 around the phosphorylation site of Tyrosine 1284.   |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000  |

## Anti-ErbB-4 / HER4 (phospho-Y1284) Antibody (A00296Y1284) Images



WesternBlot (WB) analysis of ErbB-4 / HER4 (phospho-Y1284) polyclonal antibody

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ErbB-4 / HER4 (phospho-Y1284) Antibody

For Research Use Only. Not for use in diagnostic procedures.