

## Anti-Pro-opiomelanocortin POMC Antibody

Catalog Number: A00305-1

### Overview

Product Name	Anti-Pro-opiomelanocortin POMC Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Pro-opiomelanocortin POMC Antibody catalog # A00305-1. Tested in IHC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01189

### Technical Details

Immunogen	Synthesized peptide derived from human Adrenocorticotropin(ACTH)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC-p 1:50-200, WB 1:500-2000

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.

