

Anti-PI3K p85 Alpha /p55 Gamma (phospho-Y467/Y199) PIK3R1 Antibody

Catalog Number: A00318Y467Y199


Overview

Product Name	Anti-PI3K p85 Alpha /p55 Gamma (phospho-Y467/Y199) PIK3R1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PI3K p85 Alpha /p55 Gamma (phospho-Y467/Y199) PIK3R1 Antibody catalog # A00318Y467Y199. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P27986/Q92569

Technical Details

Immunogen	Synthetic phosphopeptide derived from human PI3K p85alpha/p55gamma around the phosphorylation site of Tyrosine 467/199.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-PI3K p85 Alpha /p55 Gamma (phospho-Y467/Y199) PIK3R1 Antibody (A00318Y467Y199) Images



Western blot (WB) analysis of PI3K p85alpha/p55gamma (phospho-Y467/Y199) polyclonal antibody polyclonal antibody at 1:500 dilution Lane1:PC12 whole cell lysate Lane2:sp2/0 whole cell lysate Lane3:HEK293T whole cell lysate

1 Publications Citing This Product

1. PubMed ID: -, Lu Kong, Yongya Wu, Wangcheng Hu, Lin Liu, Yuying Xue, Geyu Liang, Mechanisms underlying reproductive toxicity induced by nickel nanoparticles identified by comprehensive gene expression analysis in GC-1 spg cells, Environmental Pollution, 2021, 116556, ISSN 0269-7

Visit bosterbio.com/anti-pi3k-p85-alpha-p55-gamma-phospho-y467-y199-pik3r1-antibody-a00318y467y199-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PI3K p85 Alpha /p55 Gamma (phospho-Y467/Y199) PIK3R1 Antibody

For Research Use Only. Not for use in diagnostic procedures.