

## Anti-PI3K p85 Alpha /p55 Gamma (Tyr607) PIK3R1 Antibody

Catalog Number: A00318Y607


### Overview

Product Name	Anti-PI3K p85 Alpha /p55 Gamma (Tyr607) PIK3R1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PI3K p85 Alpha /p55 Gamma (Tyr607) PIK3R1 Antibody catalog # A00318Y607. Tested in WB,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P27986/Q92569


### Technical Details

Immunogen	Synthetic peptide, corresponding to Human PI3K p85alpha.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:2000-1:5000 IF: 1:50-1:200

## Anti-PI3K p85 Alpha /p55 Gamma (Tyr607) PIK3R1 Antibody (A00318Y607) Images



Western blot (WB) analysis of PI3K p85alpha (Tyr607) pAb at 1:2000 dilution Lane1:C6 whole cell lysate(40ug) Lane2:CT-26 whole cell lysate(40ug) Lane3:HeLa whole cell lysate(40ug) Lane4:Jurkat whole cell lysate(40ug) Lane5:MCF-7 whole cell lysate(40ug)



Immunofluorescence analysis of BV2 cells using PI3K p85alpha/p55gamma antibody at dilution of 1:50.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PI3K p85 Alpha /p55 Gamma (Tyr607) PIK3R1 Antibody  
For Research Use Only. Not for use in diagnostic procedures.