

Anti-DNMT3B Rabbit Polyclonal Antibody

Catalog Number: A00319-1


Overview

Product Name	Anti-DNMT3B Rabbit Polyclonal Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Rabbit Polyclonal anti-DNMT3B Antibody. Catalog# A00319-1. Tested in WB. This antibody reacts with Human, Mouse, Rat.
Conjugate	Unconjugated
Application	WB
Clonality	Polyclonal
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UBC3


Technical Details

Immunogen	Recombinant protein of human DNMT3B
Isotype	IgG
Concentration	2.24 mg/ml
Purification	Purified from the immunized serum by affinity chromatography (Protein A/G)
Suggested Dilutions	WB: 1:500~2000

Anti-DNMT3B Rabbit Polyclonal Antibody (A00319-1) Images



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DNMT3B (Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DNMT3.



Western blot analysis of extracts (35ug) from different cell lines and tissues by using anti-DNMT3B rabbit polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNMT3B Rabbit Polyclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.