

Anti-Leptin Receptor/LEPR Antibody Biotin Conjugated

Catalog Number: A00350-Biotin

About LEPR

Leptin receptor, also known as LEP-R or OB-R is a protein that in humans is encoded by the LEPR gene. The protein encoded by this gene belongs to the gp130 family of cytokine receptors that are known to stimulate gene transcription via activation of cytosolic STAT proteins. This protein is a receptor for leptin (an adipocyte-specific hormone that regulates body weight), and is involved in the regulation of fat metabolism, as well as in a novel hematopoietic pathway that is required for normal lymphopoiesis. Mutations in this gene have been associated with obesity and pituitary dysfunction. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-Leptin Receptor/LEPR Antibody Biotin Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | WB, IHC, ELISA |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P48357 |

Technical Details

| | |
|---------------------|--|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human Leptin Receptor, different from the related mouse sequence by nine amino acids, and from the related rat sequence by eight amino acids. |
| Cross Reactivity | No cross-reactivity with other proteins |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Biotin |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by |

end users.
ELISA, Optimal dilutions should be determined by end users.

2 Publications Citing This Product

1. PubMed ID: -, Immunohistochemical Localization of Leptin and Leptin-Receptor Proteins in Different Tissues of Chinese Alligator, Alligator sinensis. HuaBin Zhang, YaSen Li, Tao Pan, Peng Yan, En Li, Hui Xue, XiaoBing Wu

2. PubMed ID: 26529315, The Influence of LepR Tyrosine Site Mutations on Mouse Ovary Development and Related Gene Expression Changes

Visit bosterbio.com/anti-leptin-receptor-antibody-a00350-boster.html to see all 2 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Leptin Receptor/LEPR Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.