

Anti-Cdk6 Antibody Picoband® HRP Conjugated

Catalog Number: A00358-HRP

About CDK6

Cell division protein kinase 6, also called Plstire, is an enzyme that in humans is encoded by the CDK6 gene. The protein encoded by this gene is a member of the cyclin-dependent protein kinase (CDK) family. Radiation hybrid analysis and inclusion within a mapped clone place the CDK6 gene at 7q21. Serine/threonine-protein kinase involved in the control of the cell cycle and differentiation and promotes G1/S transition. This gene also involved in initiation and maintenance of cell cycle exit during cell differentiation. It prevents cell proliferation and regulates negatively cell differentiation, but is required for the proliferation of specific cell types. In addition, CDK6 plays a role in promoting the proliferation of beta-cells in pancreatic islets of Langerhans.

Overview

Product Name	Anti-Cdk6 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q00534

Technical Details

Immunogen	E. coli-derived human Cdk6 recombinant protein (Position: M1-A326).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by

end users.
ELISA, Optimal dilutions should be determined by end users.

5 Publications Citing This Product

1. PubMed ID: 10.1186/s13046-016-0297-z, Long noncoding RNA NEAT1 promotes laryngeal squamous cell cancer through regulating miR-107/CDK6 pathway
2. PubMed ID: 34144025, Sun D,Zhao T,Long K,Wu M,Zhang Z.Triclosan down-regulates fatty acid synthase through microRNAs in HepG2 cells.Eur J Pharmacol .2021 Jun 15:174261.doi:10.1016/j.ejphar.2021.174261.Epub ahead of print.PMID:34144025.
3. PubMed ID: 27065079, TALENs-directed?knockout?of the?full-length?transcription?factor?Nrf1?that?represses?malignantbehaviour?of?human?hepatocellular?carcinoma?(HepG2)?cells

Visit bosterbio.com/anti-cdk6-picoband-trade-antibody-a00358-boster.html to see all 5 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cdk6 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.