

## Anti-Thyroglobulin TG Antibody

Catalog Number: A00359

### Overview


Product Name	Anti-Thyroglobulin TG Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Thyroglobulin TG Antibody catalog # A00359. Tested in ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01266

### Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 2511-2560
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

## Anti-Thyroglobulin TG Antibody (A00359) Images

---



Immunohistochemistry validation of TG using Anti-Thyroglobulin TG Antibody (A00359). Immunohistochemical analysis of paraffin-embedded human-thyroid

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Thyroglobulin TG Antibody

For Research Use Only. Not for use in diagnostic procedures.