

## Anti-CD47 Antibody Picoband® Cy3 Conjugated

Catalog Number: A00360-1-Cy3

### About CD47

CD47, also known as IAP or MER6, is a transmembrane protein that in humans is encoded by the CD47 gene. CD47 gene belongs to the immunoglobulin superfamily. This gene is mapped to 3q13.12. CD47 gene encodes a membrane protein, which is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to extracellular matrix. The encoded protein is also a receptor for the C-terminal cell binding domain of thrombospondin, and it may play a role in membrane transport and signal transduction. This gene has broad tissue distribution, and is reduced in expression on Rh erythrocytes.

### Overview

Product Name	Anti-CD47 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q08722

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human CD47, which shares 61.8% and 66.7% amino acid (aa) sequence identity with mouse and rat CD47, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD47 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.