

Anti-CD47 Antibody Picoband® HRP Conjugated

Catalog Number: A00360-1-HRP

About CD47

CD47, also known as IAP or MER6, is a transmembrane protein that in humans is encoded by the CD47 gene. CD47 gene belongs to the immunoglobulin superfamily. This gene is mapped to 3q13.12. CD47 gene encodes a membrane protein, which is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to extracellular matrix. The encoded protein is also a receptor for the C-terminal cell binding domain of thrombospondin, and it may play a role in membrane transport and signal transduction. This gene has broad tissue distribution, and is reduced in expression on Rh erythrocytes.

Overview

Product Name	Anti-CD47 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q08722

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human CD47, which shares 61.8% and 66.7% amino acid (aa) sequence identity with mouse and rat CD47, respectively.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.

ELISA, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD47 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.