

Anti-PDGFR-Alpha (phospho-Y1018) PDGFRA Antibody

Catalog Number: A00366Y1018


Overview

Product Name	Anti-PDGFR-Alpha (phospho-Y1018) PDGFRA Antibody
Reactive Species	Human
Description	Boster Bio Anti-PDGFR-Alpha (phospho-Y1018) PDGFRA Antibody catalog # A00366Y1018. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16234

Technical Details

Immunogen	A synthetic phosphorylated peptide around Y1018 of human PDGFRA(NP_006197.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-PDGFR-Alpha (phospho-Y1018) PDGFRA Antibody (A00366Y1018) Images



Western blot analysis of extracts of HepG2 cells, using Phospho-PDGFR alpha-Y1018 antibody at 1:1000 dilution. HepG2 cells were treated by UV for 15-30 minutes or treated by EGF for 30 minutes after serum-starvation overnight. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDGFR-Alpha (phospho-Y1018) PDGFRA Antibody

For Research Use Only. Not for use in diagnostic procedures.