

Anti-Angiopoietin-2/ANGPT2 Antibody Picoband® Cy3 Conjugated

Catalog Number: A00370-2-Cy3

About ANGPT2

ANGPT2, also known as ANG2 or Angiopoietin 2, is a protein that in humans is encoded by the ANGPT2 gene. It is mapped to 8p23.1. ANGPT2 is a naturally occurring antagonist of ANG1 that competes for binding to the TIE2 receptor and blocks ANGPT1-induced TIE2 autophosphorylation during vasculogenesis. The encoded protein disrupts the vascular remodeling ability of ANGPT1 and may induce endothelial cell apoptosis. ANGPT2 was significantly increased in plasma and alveolar edema fluid in adults with acute lung injury compared to controls or patients with hydrostatic pulmonary edema, tracheal. ANGPT2 was also significantly increased in neonates with respiratory distress syndrome who developed bronchopulmonary edema. It is also a mediator of epithelial necrosis with an important role in hyperoxic acute lung injury and pulmonary edema.

Overview

Product Name	Anti-Angiopoietin-2/ANGPT2 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O15123

Technical Details

Immunogen	E.coli-derived human Angiopoietin-2/ANGPT2 recombinant protein (Position: N357-K404).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

6 Publications Citing This Product

1. PubMed ID: 10.1620/tjem.224.221, Exercise Improves Recovery after Ischemic Brain Injury by Inducing the Expression of Angiopoietin-1 and Tie-2 in Rats
2. PubMed ID: 10.1016/j.neulet.2016.12.066, Expression of angiopoietins in central nervous system hemangioblastomas is associated with cyst formation
3. PubMed ID: 10.1007/s00198-014-2992-y, Local renin-angiotensin system is associated with bone mineral density of glucocorticoid-induced osteoporosis patients

Visit bosterbio.com/anti-angiopoietin-2-angpt2-picoband-trade-antibody-a00370-2-boster.html to see all 6 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Angiopoietin-2/ANGPT2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.