

## Anti-ATF4 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A00371-1-Fluoro647

### About ATF4

ATF4, Activating Transcription Factor 4, is also known as CREB2. ATF4 belongs to the large ATF/CREB family of transcription factors which bind DNA via their basic region and dimerize via their leucine zipper domain to form a variety of homo- and heterodimers to regulate gene transcription. It is identified that members of this family share significant sequence similarity within a leucine zipper DNA-binding motif and an adjacent basic region. The ATF4 gene is mapped to chromosome 22. Unlike CREB, which activates transcription from CRE-containing promoters, CREB2 functions as a specific repressor of CRE-dependent transcription. The transcriptional repressor activity resides within the C-terminal leucine zipper and basic domain region of the CREB2 protein.

### Overview

Product Name	Anti-ATF4 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P18848

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human ATF4, which shares 97.6% amino acid (aa) sequence identity with both mouse and rat ATF4.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 23826386, MARVELD1 Inhibits Nonsense-Mediated RNA Decay by Repressing Serine Phosphorylation of UPF1

Visit [bosterbio.com/anti-atf4-picoband-trade-antibody-a00371-1-boster.html](https://bosterbio.com/anti-atf4-picoband-trade-antibody-a00371-1-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ATF4 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.