

## Anti-Pro-epidermal growth factor EGF Antibody

Catalog Number: A00378-2

### Overview


Product Name	Anti-Pro-epidermal growth factor EGF Antibody
Reactive Species	Human
Description	Boster Bio Anti-Pro-epidermal growth factor EGF Antibody catalog # A00378-2. Tested in WB,IHC,IP applications. This antibody reacts with Human.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P01133

### Technical Details


Immunogen	recombinant protein
Cross Reactivity	Does not cross-react with eNOS or nNOS.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:100 IP: 1:10-1:50

## Anti-Pro-epidermal growth factor EGF Antibody (A00378-2) Images

---



Immunohistochemical analysis of paraffin-embedded human kidney tissue using anti-EGF antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded human skin tissue using anti-EGF antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Pro-epidermal growth factor EGF Antibody

For Research Use Only. Not for use in diagnostic procedures.