

Anti-IL-2/IL2 Antibody Picoband®

Catalog Number: A00387-1

About IL2

IL-2, Interleukin-2, formerly referred to as T-cell growth factor, is a powerfully immunoregulatory lymphokine that is produced by lectin-or antigen-activated T cells. In situ hybridization, the IL-2 gene is assigned to 4q26-q28. IL2 can act as a growth hormone for both B and T lymphocytes. It is useful in the study of the molecular nature of T-cell differentiation and, like interferons, augments natural killer cell activity.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-IL-2/IL2 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-IL-2/IL2 Antibody Picoband® catalog # A00387-1. Tested in ELISA, IHC, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl and 0.2mg Na ₂ HPO ₄ .
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P60568

Technical Details


Immunogen	E. coli-derived human IL-2 recombinant protein (Position: A21-T153). Human IL-2 shares 57.5% and 66.4% amino acid (aa) sequence identity with mouse and rat IL-2, respectively.
Cross Reactivity	No cross-reactivity with other proteins.

Isotype	Rabbit IgG
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.1-0.5ug/ml, Human Immunohistochemistry (Paraffin-embedded Section), 2-5ug/ml, Human ELISA, 0.1-0.5ug/ml, -


Anti-IL-2/IL2 Antibody Picoband® (A00387-1) Images



Western blot analysis of IL-2 using anti-IL-2 antibody (A00387-1). Electrophoresis was performed on a 12% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. Lane 1: recombinant human IL-2 protein 25 ng, Lane 2: recombinant human IL-2 protein 10 ng, Lane 3: recombinant human IL-2 protein 5 ng, Lane 4: recombinant human IL-2 protein 1 ng. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-IL-2 antigen affinity purified polyclonal antibody (A00387-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody (Catalog # BA1054) at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for IL-2 at approximately 15 kDa.



IHC analysis of IL-2 using anti-IL-2 antibody (A00387-1). IL-2 was detected in a paraffin-embedded section of human tonsil tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-IL-2 Antibody (A00387-1) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.



IHC analysis of IL-2 using anti-IL-2 antibody (A00387-1). IL-2 was detected in a paraffin-embedded section of human spleen tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-IL-2 Antibody (A00387-1) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

4 Publications Citing This Product

1. PubMed ID: 10.1007/s13273-012-0041-5, Effects of recombinant hIL-2-Luffin P1 immunotoxin combining with arotinoid ethylester on the proliferation and apoptosis of Hut-78 cells

2. PubMed ID: 26064167, Gong T, Wang Cf, Yuan Jr, Li Y, Gu Jf, Zhao Bj, Zhang L, Jia Xb, Feng L, Liu Sl. Evid Based Complement Alternat Med. 2015;2015:630760. Doi: 10.1155/2015/630760. Epub 2015 May 3. Inhibition Of Tumor Growth And Immunomodulatory Effects Of Flavonoids...

3. PubMed ID: 27284330, Treatment of radiation-induced acute intestinal injury with bone marrow-derived mesenchymal stem cells

Visit bosterbio.com/anti-il-2-picoband-trade-antibody-a00387-1-boster.html to see all 4 publications.

Anti-IL-2/IL2 Antibody

For Research Use Only. Not for use in diagnostic procedures.