

Anti-TRAF2 (S11) Antibody

Catalog Number: A00400S11


Overview

Product Name	Anti-TRAF2 (S11) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TRAF2 (S11) Antibody catalog # A00400S11. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q12933

Technical Details

Immunogen	Synthetic peptide, corresponding to the N-terminus of Human TRAF2.
Cross Reactivity	Does not cross-react with eNOS or nNOS.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 118% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-TRAF2 (S11) Antibody (A00400S11) Images



Western blot (WB) analysis of TRAF2 (S11) polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate Lane2:Raw264.7 whole cell lysate Lane3:PC12 whole cell lysate

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRAF2 (S11) Antibody

For Research Use Only. Not for use in diagnostic procedures.