

Anti-GATA-1 (phospho-S142) Antibody

Catalog Number: A00408S142

Overview

Product Name	Anti-GATA-1 (phospho-S142) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GATA-1 (phospho-S142) Antibody catalog # A00408S142. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15976

Technical Details

Immunogen	Synthetic phosphopeptide derived from human GATA1 around the phosphorylation site of Serine 142.
Cross Reactivity	Does not cross-react with eNOS or nNOS.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:50-1:200

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GATA-1 (phospho-S142) Antibody

For Research Use Only. Not for use in diagnostic procedures.