

Anti-MMP-13 Antibody

Catalog Number: A00420-1


Overview

Product Name	Anti-MMP-13 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-MMP-13 Antibody catalog # A00420-1. Tested in WB,ICC/IF,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P45452


Technical Details

Immunogen	recombinant protein
Cross Reactivity	Does not cross-react with eNOS or nNOS.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000 ICC: 1:100-1:500 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-MMP-13 Antibody (A00420-1) Images



Western blot analysis of MMP13 on MCF-7 cells lysates using anti-MMP13 antibody at 1/1,000 dilution.



ICC staining MMP13 in HeLa cells (red). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

1 Publications Citing This Product

1. PubMed ID: 32854712, Lan Q, Lu R, Chen H, Pang Y, Xiong F, Shen C, Qin Z, Zheng L, Xu G, Zhao J. MMP-13 enzyme and pH responsive theranostic nanoplatfrom for osteoarthritis. J Nanobiotechnology. 2020 Aug 27; 18(1):117. doi:10.1186/s12951-020-00666-7. PMID:32854712; PMCID:PMC7450974.

Visit bosterbio.com/anti-mmp-13-antibody-a00420-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MMP-13 Antibody

For Research Use Only. Not for use in diagnostic procedures.