

Anti-TRAF3 Antibody

Catalog Number: A00424-2


Overview

Product Name	Anti-TRAF3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TRAF3 Antibody catalog # A00424-2. Tested in WB,ELISA,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13114


Technical Details

Immunogen	KLH conjugated synthetic peptide derived from human TRAF3: 101-200/568
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ELISA: 1:5000-1:10000 IHC: 1:100-1:500 IF: 1:100-1:500

Anti-TRAF3 Antibody (A00424-2) Images



Sample:Spleen(Mouse)Lysate at 40 ug Primary: Anti-TRAF3 at 1/300 dilution



Paraformaldehyde-fixed, paraffin embedded (Mouse brain); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (TRAF3) Polyclonal Antibody, Unconjugated at 1:400 overnight at 4°C.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRAF3 Antibody

For Research Use Only. Not for use in diagnostic procedures.