

## Anti-Transcription factor 7-like 2 TCF7L2 Antibody

Catalog Number: A00431-1


### Overview

Product Name	Anti-Transcription factor 7-like 2 TCF7L2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Transcription factor 7-like 2 TCF7L2 Antibody catalog # A00431-1. Tested in WB,ICC/IF,IHC,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NQB0


### Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

## Anti-Transcription factor 7-like 2 TCF7L2 Antibody (A00431-1) Images



Western blot analysis of TCF7L2 on Jurkat cells lysates using anti-TCF7L2 antibody at 1/1,000 dilution.



ICC staining TCF7L2 in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Transcription factor 7-like 2 TCF7L2 Antibody

For Research Use Only. Not for use in diagnostic procedures.