

Anti-DNA topoisomerase 1 TOP1 Antibody

Catalog Number: A00434-2


Overview

Product Name	Anti-DNA topoisomerase 1 TOP1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-DNA topoisomerase 1 TOP1 Antibody catalog # A00434-2. Tested in WB,ICC/IF,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P11387


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 ICC: 1:50-1:500 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-DNA topoisomerase 1 TOP1 Antibody (A00434-2) Images



Western blot analysis of TOP1 on different lysates using anti-TOP1 antibody at 1/1,000 dilution. Positive control: Lane 1: HepG2 Lane 2: Jurkat Lane 3: MCF-7



ICC staining TOP1 in HeLa cells (green). The nuclear counterstain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNA topoisomerase 1 TOP1 Antibody

For Research Use Only. Not for use in diagnostic procedures.