

Anti-ROR Gamma RORC Antibody

Catalog Number: A00444-1


Overview

Product Name	Anti-ROR Gamma RORC Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ROR Gamma RORC Antibody catalog # A00444-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P51449

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 114-159 of Human RORgamma.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-ROR Gamma RORC Antibody (A00444-1) Images



Western blot (WB) analysis of RORgamma polyclonal antibody at 1:500 dilution Lane1:Hela cell lysate Lane2:sp2/0 cell lysate Lane3:PC12 cell lysate

1 Publications Citing This Product

1. PubMed ID: 28619762, Li, X.M., Chen, X., Gu, W., Guo, Y.J., Cheng, Y., Peng, J., & Guo, X.J. (2017). Impaired TNF/TNFR2 signaling enhances Th2 and Th17 polarization and aggravates allergic airway inflammation. American Journal of Physiology - Lung Cellular and Molecul...

Visit bosterbio.com/anti-ror-gamma-rorc-antibody-a00444-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ROR Gamma RORC Antibody

For Research Use Only. Not for use in diagnostic procedures.