

Anti-NPM (phospho-T199) NPM1 Antibody

Catalog Number: A00450T199


Overview

Product Name	Anti-NPM (phospho-T199) NPM1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NPM (phospho-T199) NPM1 Antibody catalog # A00450T199. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P06748

Technical Details

Immunogen	Synthetic phosphopeptide derived from human NPM around the phosphorylation site of Threonine 199.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 97% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-NPM (phospho-T199) NPM1 Antibody (A00450T199) Images



Western blot (WB) analysis of p-NPM (T199) pAb at 1:500 dilution Lane1:H1792 whole cell lysate(40ug) Lane2:A549 whole cell lysate(40ug) Lane3:SGC7901 whole cell lysate(40ug) Lane4:EC9706 whole cell lysate(40ug) Lane5:PC12 whole cell lysate(40ug) Lane6:CT26 whole cell lysate(40ug) Lane7:HEK293T whole cell lysate(40ug)

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NPM (phospho-T199) NPM1 Antibody

For Research Use Only. Not for use in diagnostic procedures.