

Anti-IL8 beta CXCR2 Antibody

Catalog Number: A00455


Overview

Product Name	Anti-IL8 beta CXCR2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-IL8 beta CXCR2 Antibody catalog # A00455. Tested in ELISA, IHC, WB applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P25025


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CXCR2. AA range:1-50
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 ELISA 1:20000 IF 1:50-200


Anti-IL8 beta CXCR2 Antibody (A00455) Images



Western blotting validation for Anti-IL8 beta CXCR2 Antibody A00455 Western Blot (WB) analysis of HuvEc cells using IL-8Rbeta polyclonal antibody. Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration



Immunohistochemistry validation of CXCR2 using Anti-IL8 beta CXCR2 Antibody (A00455). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Brain



Immunohistochemistry validation of CXCR2 using Anti-IL8 beta CXCR2 Antibody (A00455). Immunohistochemistry (IHC) analysis of paraffin-embedded Human Brain

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IL8 beta CXCR2 Antibody

For Research Use Only. Not for use in diagnostic procedures.