

Anti-ApoA LPA Monoclonal Antibody

Catalog Number: A00472


Overview

Product Name	Anti-ApoA LPA Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-ApoA LPA Monoclonal Antibody catalog # A00472. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 8F6A9; 8H5C5
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P08519


Technical Details

Immunogen	Purified recombinant fragment of ApoA (4330-4521) expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-ApoA LPA Monoclonal Antibody (A00472) Images



Western Blot (WB) analysis using ApoA Monoclonal antibody against truncated LPA-His recombinant protein (1) and truncated Trx-LPA(aa4330-4521) recombinant protein (2).



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Liver Tissues with AEC staining using ApoA Monoclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ApoA LPA Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.