

## Anti-OGT Antibody

Catalog Number: A00474


### Overview

Product Name	Anti-OGT Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-OGT Antibody catalog # A00474. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O15294

### Technical Details

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 857-1036 of human OGT (NP_858059.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2000

## Anti-OGT Antibody (A00474) Images



Western blot analysis of extracts of various cell lines, using OGT antibody at 1:1000 dilution.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-OGT Antibody

For Research Use Only. Not for use in diagnostic procedures.