

## Anti-Leptin Antibody Picoband®

Catalog Number: A00479-3

### About Lep

Leptin have been found to be obese and diabetic and to have reduced activity, metabolism and body temperature. cDNA clones encoding Leptin have been isolated from human, simian, mouse, and rat cells. The expression of Leptin mRNA has been shown to be restricted to adipose tissue. Although regulation of fat stores is deemed to be the primary function of leptin, it also plays a role in other physiological processes, as evidenced by its multiple sites of synthesis other than fat cells, and the multiple cell types beside hypothalamic cells that have leptin receptors. Many of these additional functions are yet to be defined.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Leptin Antibody Picoband®                                                                                                                                                                                                                                                                                                                                                                                                 |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                              |
| Description          | Boster Bio Anti-Leptin Antibody Picoband® catalog # A00479-3. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance. |
| Application          | WB, ELISA (Cap)                                                                                                                                                                                                                                                                                                                                                                                                                |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Formulation          | Each vial contains 4mg Trehalose, 0.9mg NaCl, 0.2mg Na2HPO4, 0.05mg NaN3.                                                                                                                                                                                                                                                                                                                                                      |
| Storage Instructions | Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                          |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Uniprot ID           | P41159                                                                                                                                                                                                                                                                                                                                                                                                                         |

### Technical Details

|                               |                                                                                           |
|-------------------------------|-------------------------------------------------------------------------------------------|
| Immunogen                     | E. coli-derived human Leptin recombinant protein (Position: V22-C167).                    |
| Recommended Detection Systems | Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western |

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
|                     | blot.                                                                     |
| Cross Reactivity    | No cross-reactivity with other proteins.                                  |
| Isotype             | Rabbit IgG                                                                |
| Form                | Lyophilized                                                               |
| Concentration       | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |
| Purification        | Immunogen affinity purified.                                              |
| Suggested Dilutions | Western blot, 0.1-0.5ug/ml<br>ELISA (Cap), 1-5ug/ml                       |

## Anti-Leptin Antibody Picoband® (A00479-3) Images

100KD —  
70KD —  
55KD —  
35KD —  
25KD —  
15KD —

Western blot analysis of Leptin using anti-Leptin antibody (A00479-3). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. Lane 1: recombinant human Leptin protein 1ng. After Electrophoresis, proteins were transferred to a Nitrocellulose membrane at 150mA for 50-90 minutes. Blocked the membrane with 5% Non-fat Milk/ TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-Leptin antigen affinity purified polyclonal antibody (Catalog # A00479-3) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:10000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for Leptin at approximately 16KD. The expected band size for Leptin is at 16KD.

## 5 Publications Citing This Product

1. PubMed ID: 10.1002/jcp.30129, miR-375 induces adipogenesis through targeting Erk1 in pancreatic duct cells under the influence of sodium palmitate
2. PubMed ID: 33107061, Gezinci-Oktayoglu S, Sancar S, Karatug-Kacar A, Bolkent S. miR-375 induces adipogenesis through targeting Erk1 in pancreatic duct cells under the influence of sodium palmitate. J Cell Physiol. 2020 Oct 26. doi: 10.1002/jcp.30129. Epub ahead of print. PMID: 3310
3. PubMed ID: 25832621, Increased osteogenesis in osteoporotic bone marrow stromal cells by overexpression of leptin

Visit [bosterbio.com/anti-leptin-picoband-trade-antibody-a00479-3-boster.html](http://bosterbio.com/anti-leptin-picoband-trade-antibody-a00479-3-boster.html) to see all 5 publications.

### Anti-Leptin Antibody

For Research Use Only. Not for use in diagnostic procedures.