

Anti-PDX1 Antibody

Catalog Number: A00491-1


Overview

Product Name	Anti-PDX1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-PDX1 Antibody catalog # A00491-1. Tested in WB,IHC,ICC/IF applications. This antibody reacts with Human.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P52945


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	IHC: 1:50-1:200 ICC/IF: 1:50-1:200 WB: 1:500

Anti-PDX1 Antibody (A00491-1) Images



Immunohistochemical analysis of paraffin-embedded human pancreas tissue using anti-PDX1 antibody. Counter stained with hematoxylin.



IF analysis of paraffin-embedded human pancreas tissue using anti-PDX1 antibody (green). The nuclear counter stain is DAPI (blue).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PDX1 Antibody

For Research Use Only. Not for use in diagnostic procedures.