

Anti-NQO1 Antibody

Catalog Number: A00494-1


Overview

Product Name	Anti-NQO1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NQO1 Antibody catalog # A00494-1. Tested in WB,ICC/IF,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15559


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:5,000 FC: 1:50-1:100 ICC: 1:50-1:200 IP: 1:50

Anti-NQO1 Antibody (A00494-1) Images



Western blot analysis of NQO1 on different lysates using anti-NQO1 antibody at 1/1,000 dilution. Positive control: Lane 1: SH-SY-5Y Lane 2: Mouse kidney Lane 3: Raji



Flow cytometric analysis of HeLa cells with NQO1 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NQO1 Antibody

For Research Use Only. Not for use in diagnostic procedures.