

Anti-Filamin 1 (R2146) FLNA Antibody

Catalog Number: A00502-2


Overview

Product Name	Anti-Filamin 1 (R2146) FLNA Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Filamin 1 (R2146) FLNA Antibody catalog # A00502-2. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P21333


Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 2100-2200 of Human Filamin 1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200


Anti-Filamin 1 (R2146) FLNA Antibody (A00502-2) Images



Western blot (WB) analysis of Filamin 1 (R2146) polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate Lane2:Raw264.7 whole cell lysate Lane3:H9C2 whole cell lysate



Immunohistochemistry (IHC) analyzes of Filamin 1 (R2146) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.



IF image of BS1128 stained 3T3-L1 cells. The cells were 4% paraformaldehyde fixed (20 min) and then incubated in 10% normal goat serum for 1h to permeabilise the cells and block non-specific protein-protein interactions. The cells were then incubated with the antibody Filamin 1 (R2146) pAb#BS1128(1:100) at 5µg/ml overnight at +4°C. The secondary antibody (Green) was Goat anti-Rabbit IgG (H+L) -FITC#BS10950 used at a 1/200 dilution for 1h. Hoechst 33342 Staining Kit#BD5011 was used to stain the cell nuclei (blue).

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Filamin 1 (R2146) FLNA Antibody

For Research Use Only. Not for use in diagnostic procedures.