

Anti-IDH2 Antibody

Catalog Number: A00510-1


Overview

Product Name	Anti-IDH2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-IDH2 Antibody catalog # A00510-1. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P48735


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-IDH2 Antibody (A00510-1) Images



Western blot analysis of IDH2 on different cell lysate using anti-IDH2 antibody at 1/1,000 dilution. Positive control □ Line 1: Mouse heart Line 2: Rat heart Line 3: 293T



Immunohistochemical analysis of paraffin-embedded human liver tissue using anti-IDH2 antibody. Counter stained with hematoxylin.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-IDH2 Antibody

For Research Use Only. Not for use in diagnostic procedures.