

Anti-PON1 Antibody

Catalog Number: A00516-2


Overview

Product Name	Anti-PON1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PON1 Antibody catalog # A00516-2. Tested in ELISA, IHC applications. This antibody reacts with Human.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P27169

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human PON1. AA range:51-100
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Suggested Dilutions	IHC 1:50-200 ELISA 1:10000-20000 IF 1:50-200

Anti-PON1 Antibody (A00516-2) Images



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PON1 Antibody

For Research Use Only. Not for use in diagnostic procedures.