

## Anti-Androgen Receptor/AR Antibody Picoband® Biotin Conjugated

Catalog Number: A00542-Biotin

### About AR

The AR (androgen receptor) gene is more than 90 kb long and codes for a protein that has 3 major functional domains: the N-terminal domain, DNA-binding domain, and androgen-binding domain. The AR gene is mapped to Xq12. The protein functions as a steroid-hormone activated transcription factor. Upon binding the hormone ligand, the receptor dissociates from accessory proteins, translocates into the nucleus, dimerizes, and then stimulates transcription of androgen responsive genes. This gene contains 2 polymorphic trinucleotide repeat segments that encode polyglutamine and polyglycine tracts in the N-terminal transactivation domain of its protein. Expansion of the polyglutamine tract causes spinal bulbar muscular atrophy (Kennedy disease). Mutations in this gene are also associated with complete androgen insensitivity (CAIS). Two alternatively spliced variants encoding distinct isoforms have been described.

### Overview

Product Name	Anti-Androgen Receptor/AR Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P10275

### Technical Details

Immunogen	E.coli-derived human Androgen Receptor/AR recombinant protein (Position: A629-Q920).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

## 2 Publications Citing This Product

---

1. PubMed ID: 10.3892/etm.2014.2150, Elevated expression levels of androgen receptors and matrix metalloproteinase-2 and -9 in 30 cases of hepatocellular carcinoma compared with adjacent tissues as predictors of cancer invasion and staging

2. PubMed ID: 10.3980/j.issn.2222-3959.2010.01.10, Effects of extract of *Buddleja officinalis* eye drops on androgen receptors of lacrimal gland cells of castrated rats with dry eye

Visit [bosterbio.com/anti-androgen-receptor-ar-picoband-trade-antibody-a00542-boster.html](http://bosterbio.com/anti-androgen-receptor-ar-picoband-trade-antibody-a00542-boster.html) to see all 2 publications.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Androgen Receptor/AR Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.