

Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody

Catalog Number: A00542Y966


Overview

Product Name	Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody catalog # A00542Y966. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P10275


Technical Details

Immunogen	Synthesized Phospho peptide derived from human JAK2/3.at AA range: T966
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 IHC 1:50-300 ELISA 1:2000-20000


Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody (A00542Y966) Images



Western blotting validation for Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody A00542Y966 Western Blot analysis of 1) mouse liver



Western blotting validation for Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody A00542Y966 Western Blot analysis of 1



Western blotting validation for Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody A00542Y966 Western Blot analysis of hela cells using primary antibody diluted at 1:1000 (4°C overnight). Secondary antibody:Goat Anti-rabbit IgG IRDye 800 (diluted at 1:5000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-JAK2/3 (Tyr966/939) AR Antibody

For Research Use Only. Not for use in diagnostic procedures.