

Anti-ZEB1 Antibody

Catalog Number: A00548-1

Overview

Product Name	Anti-ZEB1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ZEB1 Antibody catalog # A00548-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P37275

Technical Details

Immunogen	Synthesized peptide derived from part region of human protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ZEB1 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000

1 Publications Citing This Product

1. PubMed ID: 32335297, Wang Q, Shang J, Zhang Y, Zhou Y, Tang L. MiR-451a restrains the growth and metastatic phenotypes of papillary thyroid carcinoma cells via inhibiting ZEB1. *Biomed Pharmacother.* 2020 Jul;127:109901. doi:10.1016/j.biopha.2020.109901. Epub 2020 Apr 23. PMID:32335297.

Visit [bosterbio.com/anti-zeb1-antibody-a00548-1-boster.html](https://www.bosterbio.com/anti-zeb1-antibody-a00548-1-boster.html) to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ZEB1 Antibody

For Research Use Only. Not for use in diagnostic procedures.